

**Identificazione:**

1. C6:0	8. C14:1 n-5	15. C18:0	22. C20:0	29. C22:0
2. C8:0	9. C15:0	16. C18:1 <i>cis</i> n-9	23. C20:1 n-9	30. C22:1 n-9
3. C10:0	10. C15:1 n-5	17. C18:1 <i>trans</i> n-9	24. C20:2 n-6	31. C22:2 n-6
4. C11:0	11. C16:0	18. C18:2 n-6	25. C20:3 n-6	32. C24:0
5. C12:0	12. C16:1 n-7	19. C18:3 n-6	26. C20:4 n-6	33. C22:6 n-3
6. C13:0	13. C17:0	20. C18:3 n-3	27. C20:3 n-3	34. C24:1 n-9
7. C14:0	14. C17:1 n-7	21. C18:4 n-3	28. C20:5	

**Acidi Grassi Metil Esteri  
(FAMES) C6 - C24**

Colonna: **MEGA-WAX - 0.32mm, 0.25µm, 25m**

Codice Catalogo: S-WAX-032-025-25

**Condizioni:**

Iniezione: On-Column.

Detector: FID, 250°C.

Programma di temperatura: 60°C (2min), 15°C/min, 105°C,  
3°C/min, 195°C, 5°C/min,  
220°C (1min).

Carrier Gas: Idrogeno, 30 kPa.

