

Trigliceridi Burro - (C24 - C56)

Colonna: **MEGA-1, 0.32mm ID, 0.10µm, 3m + 1m Ret. Gap Integrato**

Codice Catalogo: GAP-1-032-010-3-1

Condizioni:

Iniettore: On-Column, 1.0µL.

Detector: FID, 350°C.

Programma Temperatura: 60°C, 15°C/min, 340°C.

Carrier Gas: Idrogeno, 30kPa.

Identificazione:

- | | |
|----------------|---------|
| 1. C24 | 13. C46 |
| 2. Colesterolo | 14. C48 |
| 3. C26 | 15. C50 |
| 4. C28 | 16. C52 |
| 5. C30 | 17. C54 |
| 6. C32 | 18. C56 |
| 7. C34 | |
| 8. C36 | |
| 9. C38 | |
| 10. C40 | |
| 11. C42 | |
| 12. C44 | |

