



Triglycerides - Hazelnut Oil on MEGA-LAP Column

Column: **MEGA-LAP - 0.32mm, 0.10 μ m, 25m**

Catalog Code: C-LAP-032-010-25

MEGA-LAP (Lipid Analysis Phase) is a dedicated stationary phase optimized for Lipids, Sterols and saturated and unsaturated Triglycerides separations, stable over 370°C.

Conditions:

Injection: Split 300°C, 1.0 μ L, 1:80 Split Ratio.

Detector: FID 300°C.

Oven Temp.: 120°C, 5°C/min, 330°C.

Carrier Gas: Hydrogen, 2mL/min.

Peaks:

1. PPO
2. POS
3. POO
4. PLO + OOPo
5. PLL + PoOl
6. 800
7. 000
8. OOL
9. OOL
10. LLL

Acknowledgement: Dott. C. Mariani - Stazione Sperimentale Oli e Grassi, Via Giuseppe Colombo 79 - Milano - Italy.

