



**Carvone (+/-)**  
Chiral Separation on MEGA-DEX Columns

Column: **MEGA-DEX DAC Beta - 0.25mm, 0.25 $\mu$ m, 25m**  
Catalog Code: DEX-DAC-B-025-025-25

**Conditions:**

Injection: Split 250°C, 0.5 $\mu$ L, 1:50 Split Ratio.

Sample dilution: 1mg/1mL each in n-Hexane.

Oven Temp.: 85°C, 1.5°C/min, 160°C.

Detector: FID 250°C.

Carrier Gas: Hydrogen, 75 kPa.

**Peaks:**

1. (+/-) Carvone

