

## Trigliceridi nel Burro (C24 - C56)

Colonna: **MEGA-1 - 0.32mm, 0.10µm, 3m + 1m Ret. Gap Integrato**

Codice Catalogo: GAP-1-032-010-3-1

### Condizioni:

Iniezione: On-Column.

Detector: FID, 350°C.

Programma di temperatura: 60°C, 15°C/min, 340°C.

Carrier Gas: Idrogeno, 30kPa.

### Identificazione:

- |                |         |
|----------------|---------|
| 1. C24         | 13. C46 |
| 2. Colesterolo | 14. C48 |
| 3. C26         | 15. C50 |
| 4. C28         | 16. C52 |
| 5. C30         | 17. C54 |
| 6. C32         | 18. C56 |
| 7. C34         |         |
| 8. C36         |         |
| 9. C38         |         |
| 10. C40        |         |
| 11. C42        |         |
| 12. C44        |         |

