

## NEW MEGA-MS Stationary Phases

Â

MEGA updates its GC-MS Columns line with three new MS stationary phases.

Â

All MEGA-MS phases provide minimal bleeding and a high stability at high temperatures. Discover the new MEGA-MS stationary phases. In addition to the existing MEGA-5 MS, MEGA-1 MS and MEGA-WAX MS, we add three new MS phases:

-  
MEGA-35 MS (35% Phenyl, 65% Methyl Polysiloxane - Mid Polarity)

-  
MEGA-17 MS (50% Phenyl, 50% Methyl Polysiloxane - Mid to High Polarity)

-  
MEGA-225 MS (25% Cyanopropyl, 25% Phenyl, 50% Methyl Polysiloxane - Mid To High Polarity)